

DODATNE INFORMACIJE

MODEL: AM8517MT-FUW

OSVJETLJENJE

Svjetlo / LED	375 nm UV
Broj LED žaruljica	4 UV + 4 Bijelo

OPTIKA

Povećanje	20 ~220x
-----------	----------

SENZOR

Rezolucija	8.0 Megapixel (3840x2160)
Maksimalni broj sličica u sekundi	10 fps

KOMPATIBILNOST

Sučelje	USB 2.0
Operativni sustav	Windows 10/11 macOS 10.14 or later
Podržani formati slika (Windows)	BMP, GIF, PNG, MNG, TIF, TGA, PCX, WBMP, JP2, JPC, JPG, PGX, RAS, PNM
Podržani video formati (Windows)	WMV, FLV, SWF

KUĆIŠTE

Dimenzije	130 mm (H) x 32 mm (D)
Masa	118 g

ZNAČAJKE

Mjerenje	Da
Kalibracija	Da
Microtouch senzor	Da



INFORMACIJE

Sadržaj pakiranja	Microscope, carry pouch, calibration target, user manual, Microscope, carry pouch, calibration target, user manual, N3C-C (Close Cap) N3C-EB (Extended Open Black Cap) N3C-L (Long Cap) N3C-OB (Black Open Cap) N3C-S (Sidelight Cap)
Jamstvo	2 godine
Regulatorna odobrenja	CE, FCC, ROHS